




MODEL : Linear Abraser 5750



Pneumatic Operated Life Testers

## Application

Simulates abrasion & fatigue life testing.

## TMA Q400 Spec.

Model:	DMA Q800
Manufacturer:	Taber Industries
Load Capacity :	3 N to 20 N
Operating Stroke Length :	0.5 " to 4"
Operating Speed:	2,15,25,30,40,60 cycles/min
Abradants / Wearersers :	Wide selection available (CS 10 ,H18 , wearersers)

## Fatigue Life Testers

Custom-made Fixtures
Max Pressure : 110psi
Programmable Test counter
Pneumatic cylinders of varying stroke lengths / diameters
Flow-control valves for actuation speed control